

Elta SilverGel™

Advanced Ionic-Silver Wound Gel

Product Information

- Crystal clear, broad spectrum antimicrobial wound gel
- Amorphous gel allows full antimicrobial contact with entire wound bed
- Sustained & effective silver ion release for 3* days+
- Remains clear while in use – will not discolor or stain tissue
- Non-cytotoxic
- Ideal viscosity – very easy to apply & remains in wound bed
- Use on:
 - pressure ulcers (Stage I-IV)
 - partial and full thickness wounds
 - leg ulcers
 - diabetic foot ulcers
 - arterial ulcers
 - chronic wounds
 - grafted wounds and donor sites
 - surgical wounds
 - lacerations and abrasions
 - 1st and 2nd degree burns with light to moderate exudate

Why a wound gel?

When used appropriately, Elta Silver Wound Gel forms a coating within the wound bed. Research has shown that antimicrobial dressings containing silver must come fully into contact with the wound bed in order to be most effective:

“One additional factor that impinges on antiseptic efficacy and is not related to the form of silver used or the dosage but should not be overlooked is the ability of the carrier dressing to conform to the wound bed. High conformability helps ensure that areas of non-contact between the dressing and the wound bed are minimized thus reducing the formation of voids (dead space) where bacteria may flourish.”*

**Exploring the Effects of Silver in Wound Management - What is Optimal?*
Richard White, PhD; Keith Cutting, MN, RN, Dip N, Cert Ed
Wounds. 2006;18(11):307-314. ©2006 Health Management Publications, Inc.

Swiss-American

2055 Luna Road, #126, Carrollton, TX 75006

(800) 633-8872

www.elta.net fax: (972) 385-7930

elta®

SKIN AND WOUND CARE



General Description and Uses:

A clear, greaseless, antimicrobial hydrogel used create or maintain a moist wound environment. Used for pressure and arterial ulcers, burns, cuts & other skin irritations.

Elta SilverGel is a broad spectrum antimicrobial – 100% effective at contact area throughout a 1-3 day period of time.

Directions:

Cleanse wound with cleansing liquid, such as Elta Wound Cleanser. Evenly apply Elta SilverGel approximately 1/8 to 1/4" thick to wound area. Cover wound with appropriate dressing. Repeat as necessary to keep the wound moist.

Package size:

1 oz, 1.5 oz, 8oz and 16 oz

Applications per Package:

Varies based on size of wound. Sealable container allows multiple applications

Availability:

Contact Swiss-American or your local drug or medical supply distributor

Do not use on:

Wounds with high amount of exudate.
Patients with known sensitivity to silver.